

Product datasheet for **MP211040**

Plau Mouse qPCR Primer Pair (NM_008873)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	18792
Forward Sequence:	AGAAGCGACCCTGGTGCTATGT
Reverse Sequence:	CCACACTGGAAGCCTTGTTGGT
ACCN:	<u>BC120709</u> , <u>BC120713</u> , <u>NM_008873</u> , <u>NM_008873.1</u> , <u>NM_008873.2</u> , <u>NM_008873.3</u> , <u>BC145037</u> , <u>BY320655</u>
UniProt ID:	<u>P06869</u>
Synonyms:	u-PA; uPA
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

