

Product datasheet for **MP211050**

Pafah2 Mouse qPCR Primer Pair (NM_133880)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	100163
Forward Sequence:	CCAAGCATCAGAGAAGTGCGAG
Reverse Sequence:	CCAGCTAACAGGCAAGCGACAA
ACCN:	<u>BC021890</u> , <u>BC025495</u> , <u>NM_133880</u> , <u>NM_133880.1</u> , <u>NM_133880.2</u> , <u>NM_133880.3</u> , <u>BC025938</u> , <u>BQ179290</u> , <u>BU757520</u> , <u>BY032695</u> , <u>BY233826</u> , <u>BY782149</u>
UniProt ID:	<u>Q8VDG7</u>
Synonyms:	2310074E22Rik; AI507170
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)