

Product datasheet for **MP211060**

Psd Mouse qPCR Primer Pair (NM_028627)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	73728
Forward Sequence:	GAGGTGGAAGAGTCAGATGCTG
Reverse Sequence:	GAGTCTGAGTCGGATACCACT
ACCN:	<u>NM_028627</u> , <u>NM_028627.1</u> , <u>NM_028627.2</u> , <u>BC145352</u> , <u>BC043135</u> , <u>BC057353</u> , <u>BC058352</u> , <u>BC069847</u> , <u>BC138652</u>
UniProt ID:	<u>Q5DTT2</u>
Synonyms:	1110007H17Rik; A930015K15; Efa6; Efa6a; mKIAA2011; Psdl
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

