

Product datasheet for **MP211065**

Plekha4 Mouse qPCR Primer Pair (NM_148927)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	69217
Forward Sequence:	CTCTGGTTCGTCCACCTATCTC
Reverse Sequence:	GCTCTGATGCAAAGCTGAGTCC
ACCN:	<u>BC021387</u> , <u>BC024961</u> , <u>NM_148927</u> , <u>NM_148927.1</u> , <u>NM_148927.2</u> , <u>BC030731</u> , <u>BC034903</u> , <u>BE652135</u> , <u>BG081030</u> , <u>BY235943</u>
Synonyms:	2410005C22Rik; C88169; Pepp1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

