

Product datasheet for **MP211067**

Plekha5 Mouse qPCR Primer Pair (NM_144920)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	109135
Forward Sequence:	GACGTACAGGTCAGAAGTGACC
Reverse Sequence:	GCTGGCATTGACCTTCTGTGAG
ACCN:	<u>BC014853</u> , <u>BC021505</u> , <u>BC033481</u> , <u>NM_144920</u> , <u>NM_144920.1</u> , <u>NM_144920.2</u> , <u>NM_144920.3</u> , <u>BC014853.1</u> , <u>BC010246</u> , <u>BC025596</u> , <u>BC031390</u>
Synonyms:	2810431N21Rik; AI428202; AK129423; AMH-Cre; Ayu21-9; Gt(pU21)9Imeg; Pepp2; Tg(AMH-cre)1Flor
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)