

Product datasheet for **MP211078**

Plekhg5 Mouse qPCR Primer Pair (NM_001004156)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	269608
Forward Sequence:	ACCAAGGAAGCCATCGTCACCA
Reverse Sequence:	TCCACAACCTCGTAGGCATCGA
ACCN:	<u>NM_001004156</u> , <u>NM_001004156.1</u> , <u>NM_001004156.2</u> , <u>BC023181</u> , <u>BC064091</u> , <u>BC148338</u> , <u>BC152981</u> , <u>BY001826</u> , <u>BY151802</u>
UniProt ID:	<u>Q66T02</u>
Synonyms:	BC023181; mKIAA0720; RP23-284C20.1; Syx
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)