

Product datasheet for **MP211088**

Phip Mouse qPCR Primer Pair (NM_001081216)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	83946
Forward Sequence:	GCAGGTGTTAGTAATGCCAGCC
Reverse Sequence:	CACTTCTCGGAGCATTGCTGTG
ACCN:	<u>NM_001081216</u> , <u>NM_001081216.1</u> , <u>BC009146</u> , <u>BC022102</u> , <u>BC037036</u> , <u>BC037950</u> , <u>BC039765</u> , <u>BC049950</u> , <u>BC068022</u> , <u>BC108423</u> , <u>BC160290</u>
Synonyms:	2810004D21Rik; 4632404O06Rik; Ndrp; Wdr11
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

