

Product datasheet for **MP211091**

Cyth4 Mouse qPCR Primer Pair (NM_028195)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	72318
Forward Sequence:	GCCTTTGTGGACTGCCATGAGT
Reverse Sequence:	GCCTCCATCATGCGGTCAATCT
ACCN:	<u>NM_028195</u> , <u>NM_028195.1</u> , <u>NM_028195.2</u> , <u>NM_028195.3</u> , <u>BC139054</u> , <u>BC018505</u> , <u>BC039640</u> , <u>BC139056</u> , <u>BC145528</u> , <u>BY161747</u>
UniProt ID:	<u>Q80YW0</u>
Synonyms:	2510004M07Rik; 5830469K17Rik; AI467541; mFLJ00017; Pscd4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)