

Product datasheet for **MP211110**

Plxna4 Mouse qPCR Primer Pair (NM_175750)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	243743
Forward Sequence:	GCGGATTGAACCAGAATGGAGC
Reverse Sequence:	TTCCTCCATGCTTGGCACGGAT
ACCN:	<u>NM_175750, NM_175750.1, NM_175750.2, NM_175750.3, BC141298, BF456700</u>
UniProt ID:	<u>Q80UG2</u>
Synonyms:	9330117B14; mKIAA1550; Plxa4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)