

Product datasheet for **MP211116**

Plxdc2 Mouse qPCR Primer Pair (NM_026162)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	67448
Forward Sequence:	GTCTCCAGTTCAATGGTTGTGGC
Reverse Sequence:	GATGGCGATCAAATCCACTGGAG
ACCN:	BC057881 , NM_026162 , NM_026162.1 , NM_026162.2 , NM_026162.3 , NM_026162.4 , NM_026162.5 , NM_026162.6
UniProt ID:	Q9DC11
Synonyms:	1200007L24Rik; 5430431D22Rik; AI646728; AU022916; AV231634; Tem7r
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)