

Product datasheet for **MP211126**

Pvr Mouse qPCR Primer Pair (NM_027514)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	52118
Forward Sequence:	GAGGCAGTAGAAGCACCAATGC
Reverse Sequence:	GGTGACCATTGGCAGAGATGCA
ACCN:	<u>BC013673</u> , <u>BC032283</u> , <u>NM_027514</u> , <u>NM_027514.1</u> , <u>NM_027514.2</u> , <u>BC138219</u> , <u>BC138220</u> , <u>BY276828</u>
Synonyms:	3830421F03Rik; CD155; D7Ert458e; HVED; mE4; necl-5; PVS; Taa1; Tage4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)