

Product datasheet for **MP211129**

Nectin3 Mouse qPCR Primer Pair (NM_021495)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	58998
Forward Sequence:	GCAAAGCCGTTACATTCCCACTT
Reverse Sequence:	CACAAACGGCTGCTACTGTCTC
ACCN:	<u>BC125588</u> , <u>BC132187</u> , <u>NM_021495</u> , <u>NM_021495.1</u> , <u>NM_021495.2</u> , <u>NM_021495.3</u> , <u>NM_021495.4</u> , <u>BC096530</u> , <u>BQ174425</u> , <u>BU701648</u> , <u>BY713810</u>
UniProt ID:	<u>Q9JLB9</u>
Synonyms:	2610301B19Rik; 3000002N23Rik; 4921513D19Rik; AA407785; AU016832; AW538082; CD113; PvrI3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)