

## Product datasheet for **MP211131**

### Plk1 Mouse qPCR Primer Pair (NM\_011121)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	18817
Forward Sequence:	CCATCTTCTGGGTCAGCAAGTG
Reverse Sequence:	CCGTCATTGTAGAGAATCAGGCG
ACCN:	<u><a href="#">BC006880</a>, <a href="#">NM_011121</a>, <a href="#">NM_011121.1</a>, <a href="#">NM_011121.2</a>, <a href="#">NM_011121.3</a>, <a href="#">NM_011121.4</a>, <a href="#">BG065076</a></u>
UniProt ID:	<u><a href="#">Q07832</a></u>
Synonyms:	P; Plk; STPK; STPK13
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)