

Product datasheet for **MP211140**

Parp12 Mouse qPCR Primer Pair (NM_172893)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	243771
Forward Sequence:	AAGTTCTGACCTGGTGAGCAGG
Reverse Sequence:	TGGTGACACAGGACCTGAACTG
ACCN:	<u>BC120733</u> , <u>NM_172893</u> , <u>NM_172893.1</u> , <u>NM_172893.2</u> , <u>NM_172893.3</u> , <u>BC024579</u> , <u>BC028906</u> , <u>BC048927</u> , <u>BC137645</u>
UniProt ID:	<u>Q8BZ20</u>
Synonyms:	9930021O16; AA409132; AA536654; ARTD12; PARP-12; Zc3hdc1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

