

## Product datasheet for **MP211141**

### Parp14 Mouse qPCR Primer Pair (NM\_001039530)

#### Product data:

|                          |                                                                                                                                                                                                                                                                                                                                               |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type:</b>     | qPCR Primer Pairs                                                                                                                                                                                                                                                                                                                             |
| <b>Gene ID:</b>          | 547253                                                                                                                                                                                                                                                                                                                                        |
| <b>Forward Sequence:</b> | GGTCTCTGTTGATGTGGGCACA                                                                                                                                                                                                                                                                                                                        |
| <b>Reverse Sequence:</b> | TTCTCGGTGGTCTTCAGGCAGT                                                                                                                                                                                                                                                                                                                        |
| <b>ACCN:</b>             | <u>NM_001039530</u> , <u>NM_001039530.1</u> , <u>NM_001039530.2</u> , <u>NM_001039530.3</u> , <u>BC138874</u> , <u>BC021340</u> , <u>BC047386</u> , <u>BC145452</u>                                                                                                                                                                           |
| <b>UniProt ID:</b>       | <u>Q2EMV9</u>                                                                                                                                                                                                                                                                                                                                 |
| <b>Synonyms:</b>         | 1600029O10Rik; ARTD8; BC021340; CoaSt6; mKIAA1268                                                                                                                                                                                                                                                                                             |
| <b>Components:</b>       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions).<br>Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.                                                                                                                                |
| <b>Quality Control:</b>  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| <b>Storage:</b>          | Store at -20°C.                                                                                                                                                                                                                                                                                                                               |
| <b>Stability:</b>        | The primer mix is stable for one year from date of shipping.                                                                                                                                                                                                                                                                                  |

