

Product datasheet for **MP211145**

Parp6 Mouse qPCR Primer Pair (NM_029922)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	67287
Forward Sequence:	GAAGAAGTGGCTACAGGAGCAG
Reverse Sequence:	TCCACCACAGAGGGATAAGGCT
ACCN:	BC031387 , BC062096 , NM_029922 , NM_029922.2 , NM_029922.3 , BC062096.1 , BC030184 , BC052064 , BC058266 , BQ442721
UniProt ID:	Q6P6P7
Synonyms:	1700119G14Rik; 2310028P13Rik; 3110038K10Rik; ARTD17; C030013N01Rik; PARP-6
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)