

Product datasheet for **MP211154**

Pap0lg Mouse qPCR Primer Pair (NM_172555)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	216578
Forward Sequence:	TGGTGTCTCTGGGCAATGCTA
Reverse Sequence:	GCTTCAACAGCACAGGATTTGGC
ACCN:	BC059271 , NM_172555 , NM_172555.1 , NM_172555.2 , BC019635 , BC023030 , BC037439 , BC049163
UniProt ID:	Q6PCL9
Synonyms:	9630006B20Rik
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)