

Product datasheet for **MP211166**

Pbrm1 Mouse qPCR Primer Pair (NM_001081251)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	66923
Forward Sequence:	AACGGGCAAGAAGGATGAACCG
Reverse Sequence:	GCTGAGTGCTGGAGAAAGAAGG
ACCN:	<u>NM_001081251</u> , <u>NM_001081251.1</u> , <u>BC023452</u> , <u>BC024438</u> , <u>BC025605</u> , <u>BC029037</u> , <u>BC037018</u> , <u>BC055456</u> , <u>BC055708</u> , <u>BC172609</u> , <u>BC172736</u>
Synonyms:	2610016F04Rik; AI507524; BAF180; Pb1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

