

Product datasheet for **MP211169**

Pcgf3 Mouse qPCR Primer Pair (NM_172716)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	69587
Forward Sequence:	GGCTACCTCATTGACGCAACCA
Reverse Sequence:	GGTGGCTCTGATGGATCACGAT
ACCN:	<u>BC049266, NM_172716, NM_172716.1, NM_172716.2, NM_172716.3, NM_172716.4</u>
UniProt ID:	<u>Q8BTQ0</u>
Synonyms:	2310035N15Rik; AI662857; D630042K08Rik; DONG1; E430039C14; Rnf3; RNF3A
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)