

Product datasheet for **MP211205**

Poldip3 Mouse qPCR Primer Pair (NM_178627)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	73826
Forward Sequence:	CGAACCAAAGCCTTGACCAGCA
Reverse Sequence:	TCACGGTCATCTTGGTGCCTTC
ACCN:	<u>BC037652</u> , <u>BC046505</u> , <u>BC058225</u> , <u>NM_178627</u> , <u>NM_178627.1</u> , <u>NM_178627.2</u> , <u>NM_178627.3</u> , <u>BB711286</u> , <u>BC029751</u> , <u>BY298213</u> , <u>NM_178627.4</u>
UniProt ID:	<u>Q8BG81</u>
Synonyms:	I110008P04Rik; AA408269; AL022852; C77954; mKIAA1649; PDIP46; SKAR
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

