

Product datasheet for **MP211215**

Polr2h Mouse qPCR Primer Pair (NM_145632)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	245841
Forward Sequence:	CCTTGGATGATGGTGAATACAACC
Reverse Sequence:	GCTGCTTCAGTGGAGGTCTCAT
ACCN:	BC002306 , BC092087 , BC125612 , BC125638 , NM_145632 , NM_145632.1 , NM_145632.2 , BC092087.1
UniProt ID:	Q923G2
Synonyms:	MGC8248; MGC103049
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)