

Product datasheet for **MP211229**

Cavin1 Mouse qPCR Primer Pair (NM_008986)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	19285
Forward Sequence:	GCGCCGCAACTTCAAAGTCATG
Reverse Sequence:	TCACCCTCCTTCTCAGGCAGTG
ACCN:	BC110698 , BC116701 , BC117527 , NM_008986 , NM_008986.1 , NM_008986.2 , BC065792 , BC099493 , BY132697
UniProt ID:	O54724
Synonyms:	2310075E07Rik; AW546441; Cav-p60; Cavin
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

