

## Product datasheet for **MP211232**

### **Ptbp1 Mouse qPCR Primer Pair (NM\_001077363)**

#### **Product data:**

<b>Product Type:</b>	qPCR Primer Pairs
<b>Gene ID:</b>	19205
<b>Forward Sequence:</b>	CACCGCTTCAAGAAACCAGGCT
<b>Reverse Sequence:</b>	GTTGCTGGAGAAGAGGCTCTTG
<b>ACCN:</b>	<a href="#">BC007472</a> , <a href="#">BC028848</a> , <a href="#">BC066210</a> , <a href="#">BC086489</a> , <a href="#">NM_001077363</a> , <a href="#">NM_001077363.1</a> , <a href="#">NM_001077363.2</a> , <a href="#">BC006069</a> , <a href="#">BC006666</a> , <a href="#">BC008531</a> , <a href="#">BY167466</a>
<b>Synonyms:</b>	AA407203; AL033359; HNRPI; pPTB; Ptb; PTB-1; PTB2; PTB3; PTB4
<b>Components:</b>	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
<b>Quality Control:</b>	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix is stable for one year from date of shipping.



[View online »](#)