

Product datasheet for **MP211312**

Kcne2 Mouse qPCR Primer Pair (NM_134110)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	246133
Forward Sequence:	CGTCATCCTGTACCTCATGGTG
Reverse Sequence:	CCTGCCAATCTTCCACGATGTAC
ACCN:	<u>BC022699</u> , <u>NM_134110</u> , <u>NM_134110.1</u> , <u>NM_134110.2</u> , <u>NM_134110.3</u>
UniProt ID:	<u>Q9D808</u>
Synonyms:	2200002116Rik; AW048273; MiRP1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

