

Product datasheet for **MP211345**

Kcnq5 Mouse qPCR Primer Pair (NM_023872)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	226922
Forward Sequence:	GCATCAGTAGGTGACAGGAGATC
Reverse Sequence:	TCGGCGTCTATCACTGGCTTTG
ACCN:	<u>NM_023872</u> , <u>NM_023872.1</u> , <u>NM_023872.2</u> , <u>NM_023872.3</u> , <u>BC156308</u> , <u>BC157041</u> , <u>BE648243</u>
UniProt ID:	<u>Q9JK45</u>
Synonyms:	7730402H11; 9230107O05Rik; AA589396; D1Mgi1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)