

Product datasheet for **MP211352**

Pou3f2 Mouse qPCR Primer Pair (NM_008899)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	18992
Forward Sequence:	AACGGCGGTTTGCTCTATTCGC
Reverse Sequence:	TCTGCATGGTGTGGCTCATCGT
ACCN:	<u>NM_008899</u> , <u>NM_008899.1</u> , <u>NM_008899.2</u> , <u>BC172622</u> , <u>BC172726</u>
UniProt ID:	<u>P31360</u>
Synonyms:	9430075J19Rik; A230098E07Rik; Brn-2; Brn2; oct-7; OTF-7; Otf7
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)