

Product datasheet for **MP211363**

Slc66a2 Mouse qPCR Primer Pair (NM_025861)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	66943
Forward Sequence:	TGTCCATCGACTCGGCTCTGTT
Reverse Sequence:	CACTTGTCCACATGAGCACCATC
ACCN:	BC008231 , BC029179 , BC043686 , NM_025861 , NM_025861.1 , NM_025861.2 , NM_025861.3 , BB728808 , BB842599 , BB843279 , BY263686
UniProt ID:	Q80XM9
Synonyms:	2310009N05Rik; 4933425L21Rik; 5730564E11Rik; C78974
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

