

Product datasheet for **MP211370**

Prdm2 Mouse qPCR Primer Pair (NM_001081355)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	110593
Forward Sequence:	AGCTGCGATTGAGGAAGAGCGA
Reverse Sequence:	CCTGGGTTCTGACGCCTTTGTT
ACCN:	<u>NM_001081355</u> , <u>NM_001081355.1</u> , <u>NM_001081355.2</u> , <u>NM_001081355.3</u> , <u>BB556664</u> , <u>BC043456</u> , <u>BC052515</u> , <u>BC151062</u> , <u>BC156279</u> , <u>BU756465</u> , <u>BY232351</u>
Synonyms:	4833427P12Rik; E330024L24Rik; Gm1033; Gm19732; KMT8; Riz; Rizl; Znfprlcl
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

