

Product datasheet for **MP211372**

Prdm5 Mouse qPCR Primer Pair (NM_027547)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	70779
Forward Sequence:	GATTGGCTACCTGGACAGTGATG
Reverse Sequence:	CACACTGTGGACAGGCAAAGTC
ACCN:	<u>BC089598</u> , <u>NM_027547</u> , <u>NM_027547.1</u> , <u>NM_027547.2</u> , <u>BC138901</u> , <u>BC138902</u> , <u>BG070775</u> , <u>BX511446</u> , <u>BY014014</u>
UniProt ID:	<u>Q9CXE0</u>
Synonyms:	4432417F03Rik; 6530401I24Rik; AII9729I; EI301I2L17Rik; PFM2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

