

## Product datasheet for **MP211669**

### Prnp Mouse qPCR Primer Pair (NM\_011170)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	19122
Forward Sequence:	CAGCAACCAGAACAACCTTCGTGC
Reverse Sequence:	CGCTCCATCATCTTCACATCGG
ACCN:	<u><a href="#">BC006703</a></u> , <u><a href="#">NM_011170</a></u> , <u><a href="#">NM_011170.1</a></u> , <u><a href="#">NM_011170.2</a></u> , <u><a href="#">NM_011170.3</a></u> , <u><a href="#">BC032085</a></u> , <u><a href="#">BQ044630</a></u>
UniProt ID:	<u><a href="#">P04925</a></u>
Synonyms:	AA960666; AI325101; CD230; Prn-i; Prn-p; PrP; prP27-30; prP33-35C; PrP <sup>PrPC</sup> ; PrP <sup>Sc</sup> ; Sinc
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)