

Product datasheet for **MP211697**

Paqr3 Mouse qPCR Primer Pair (NM_198422)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	231474
Forward Sequence:	GGTTTGGCTCAATGGAGGAGTC
Reverse Sequence:	GGCTTGATCCGAGGTAGTTCAG
ACCN:	<u>BC096380</u> , <u>BC120853</u> , <u>BC120855</u> , <u>NM_198422</u> , <u>NM_198422.1</u> , <u>NM_198422.2</u> , <u>BC022731</u> , <u>BC039974</u> , <u>BC050248</u> , <u>BC051173</u> , <u>NM_198422.3</u>
UniProt ID:	<u>Q6TCG8</u>
Synonyms:	6330415A20Rik; AY424292; RKTG
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

