

Product datasheet for **MP211705**

Pdcd1lg2 Mouse qPCR Primer Pair (NM_021396)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	58205
Forward Sequence:	CTGGGACTACAAGTACCTGACG
Reverse Sequence:	CTCTAGCCTGGCAGGTAAGCTG
ACCN:	<u>BC104138</u> , <u>BC104139</u> , <u>NM_021396</u> , <u>NM_021396.1</u> , <u>NM_021396.2</u> , <u>BY749628</u>
UniProt ID:	<u>Q9WUL5</u>
Synonyms:	B7-DC; Btdc; F730015O22Rik; PD-L2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)