

Product datasheet for **MP211713**

Pdcd11 Mouse qPCR Primer Pair (NM_011053)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	18572
Forward Sequence:	AGGAGTATGGCGTGTTCGTCCA
Reverse Sequence:	GGTCACCTTAGCAACGACTGTC
ACCN:	<u>BC070468</u> , <u>NM_011053</u> , <u>NM_011053.1</u> , <u>NM_011053.2</u> , <u>BC010844</u> , <u>BC038503</u> , <u>BC051231</u> , <u>BC055276</u> , <u>BC094244</u> , <u>BY729453</u>
UniProt ID:	<u>Q6NS46</u>
Synonyms:	1110021I22Rik; ALG-4; mKIAA0185; Pdcd7
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

