

Product datasheet for **MP211725**

Prl2b1 Mouse qPCR Primer Pair (NM_025532)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	66392
Forward Sequence:	CACCATCCAGGACTCCTTAGAG
Reverse Sequence:	GTGACAAGTCCCAGGGATTCTG
ACCN:	BC100329 , BC111112 , NM_025532 , NM_025532.1 , NM_025532.2 , BC138199
UniProt ID:	Q9DAZ2
Synonyms:	2310047B08Rik; Ghd12; mPLP-K; PLP-K; Plpk; Prlpk
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

