

Product datasheet for **MP211726**

Prl2c2 Mouse qPCR Primer Pair (NM_031191)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	18811
Forward Sequence:	TTCCTTCCAACCTCCAGAAAACAAG
Reverse Sequence:	CTAGATCGTCCAGAGGGCTTTC
ACCN:	<u>BC056203</u> , <u>BC064771</u> , <u>BC080826</u> , <u>BC100300</u> , <u>NM_031191</u> , <u>NM_031191.1</u> , <u>BC145439</u> , <u>NM_031191.2</u>
UniProt ID:	<u>P04095</u>
Synonyms:	Ghd2; MRP-1; Plf; PLF-1; Plf1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

