

Product datasheet for **MP211758**

Prr11 Mouse qPCR Primer Pair (NM_175563)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	270906
Forward Sequence:	ACCAGAGGAACGCACTCCACTA
Reverse Sequence:	CTACCTGGACTTTTCTCAGCAGC
ACCN:	BC023918 , NM_175563 , NM_175563.1 , NM_175563.2 , NM_175563.3 , NM_175563.4 , NM_175563.5
UniProt ID:	Q8BHE0
Synonyms:	B930067F20Rik
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

