

Product datasheet for **MP211773**

Psrc1 Mouse qPCR Primer Pair (NM_019976)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	56742
Forward Sequence:	GGAGTCTACCACCTGAAAAGG
Reverse Sequence:	GCACTGGTGGTTGGAAGGTTTG
ACCN:	<u>BC008994</u> , <u>BC019165</u> , <u>NM_019976</u> , <u>NM_019976.1</u> , <u>NM_019976.2</u> , <u>NM_019976.3</u> , <u>BY080522</u> , <u>BY231608</u>
UniProt ID:	<u>Q9D0P7</u>
Synonyms:	5430413I02Rik; AL033312; Dda3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

