

Product datasheet for **MP211829**

Ptov1 Mouse qPCR Primer Pair (NM_133949)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	84113
Forward Sequence:	TTCACGCCTGGTACAGTCCAC
Reverse Sequence:	ACCTCACACGATGCCTTGTAGG
ACCN:	BC008543 , BC024632 , NM_133949 , NM_133949.1 , BC038159 , BC062108 , NM_133949.2
UniProt ID:	Q91VU8
Synonyms:	1110030G05Rik; AU041779; Gcap3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

