

## Product datasheet for **MP211867**

### **Psm2 Mouse qPCR Primer Pair (NM\_134101)**

#### **Product data:**

<b>Product Type:</b>	qPCR Primer Pairs
<b>Gene ID:</b>	21762
<b>Forward Sequence:</b>	CAGAGCCATTCCGCAGTTTGC
<b>Reverse Sequence:</b>	GTGCTCACTACAGATGTGGAGG
<b>ACCN:</b>	<a href="#">BC021508</a> , <a href="#">BC023049</a> , <a href="#">BC024830</a> , <a href="#">BC027388</a> , <a href="#">BC027767</a> , <a href="#">BC027822</a> , <a href="#">BC028851</a> , <a href="#">BC031128</a> , <a href="#">NM_134101</a> , <a href="#">NM_134101.1</a> , <a href="#">NM_134101.2</a> , <a href="#">BC005592</a> , <a href="#">BC024058</a> , <a href="#">BC033325</a> , <a href="#">BY307723</a>
<b>UniProt ID:</b>	<a href="#">Q8VDM4</a>
<b>Synonyms:</b>	9430095H01Rik; AA407121; TEG-190; Tex190
<b>Components:</b>	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
<b>Quality Control:</b>	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix is stable for one year from date of shipping.



[View online »](#)