

Product datasheet for **MP211873**

Psmid8 Mouse qPCR Primer Pair (NM_026545)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	57296
Forward Sequence:	GCCAAGGACATCCAGACCAATG
Reverse Sequence:	AGAAGGTGTAGCTTTCGGCAGG
ACCN:	<u>BC004075, BC005717, NM_026545, NM_026545.1, NM_026545.2, NM_026545.3, BC004075.1</u>
UniProt ID:	<u>Q9CX56</u>
Synonyms:	6720456J22Rik; AA407360; AL033291; AL033322; AL033323; C76433
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)