

Product datasheet for **MP211904**

Pin4 Mouse qPCR Primer Pair (NM_027181)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	69713
Forward Sequence:	GGGATGAGATTCACTGAAGTGGC
Reverse Sequence:	TCCATCCCACTTACAGGCAAGG
ACCN:	<u>BC119355</u> , <u>BC119357</u> , <u>NM_027181</u> , <u>NM_027181.1</u> , <u>NM_027181.2</u>
UniProt ID:	<u>Q9CWW6</u>
Synonyms:	2410002I22Rik; EPVH; Par14
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)