

Product datasheet for **MP211923**

Pias2 Mouse qPCR Primer Pair (NM_008602)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	17344
Forward Sequence:	CTTCGAGTGCCTTGATGTGCC
Reverse Sequence:	AGGACAGATCCAGGTGGGCTTT
ACCN:	<u>BC005596</u> , <u>BC034711</u> , <u>NM_008602</u> , <u>NM_008602.1</u> , <u>NM_008602.2</u> , <u>NM_008602.3</u> , <u>NM_008602.4</u> , <u>BU704371</u>
UniProt ID:	<u>Q8C5D8</u>
Synonyms:	6330408K17Rik; AI462206; ARIP3; AU018068; Dib; Miz1; PIASxalpha; PIASxb; PIASxbeta; SIZ2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)