

Product datasheet for **MP211933**

Prkcb Mouse qPCR Primer Pair (NM_008855)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	18751
Forward Sequence:	CCAAGATGACGATGTGGAGTGC
Reverse Sequence:	CTCCATCACAAAGTACAGGCGG
ACCN:	BC127083 , NM_008855 , NM_008855.1 , NM_008855.2 , BC038148 , BC048553 , BU701138 , BY729820
UniProt ID:	P68404
Synonyms:	PKC-B; PKC-Beta; Pkcb; Prkcb1; Prkcb2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)