

Product datasheet for **MP211952**

Prkaa2 Mouse qPCR Primer Pair (NM_178143)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	108079
Forward Sequence:	CTGAAGCCAGAGAATGTGCTGC
Reverse Sequence:	GAGATGACCTCAGGTGCTGCAT
ACCN:	<u>NM_178143</u> , <u>NM_178143.1</u> , <u>NM_178143.2</u> , <u>BC138566</u> , <u>BC037512</u> , <u>BC138565</u> , <u>BY006553</u> , <u>BY065557</u> , <u>BY718361</u>
UniProt ID:	<u>Q8BRK8</u>
Synonyms:	2310008I11Rik; A830082D05; AMPKalpha2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)