

Product datasheet for **MP211953**

Prkab1 Mouse qPCR Primer Pair (NM_031869)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	19079
Forward Sequence:	CCAAAAGTGCTCCGATGTGTCTG
Reverse Sequence:	GGGCTTTGAACCTCTCTTCTGG
ACCN:	<u>BC016398</u> , <u>NM_031869</u> , <u>NM_031869.1</u> , <u>NM_031869.2</u> , <u>BC027650</u>
UniProt ID:	<u>Q9R078</u>
Synonyms:	1300015D22Rik; AU021155; E430008F22
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)