

Product datasheet for **MP211958**

Prkar1b Mouse qPCR Primer Pair (NM_008923)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	19085
Forward Sequence:	ATGGCATCCAGCAGGTGCTCAA
Reverse Sequence:	CACTGGGAGTTTGACTTCTGCC
ACCN:	BC011424 , NM_008923 , NM_008923.1 , NM_008923.2 , NM_008923.3 , BB429442 , BU560122 , BU844021 , BY237508
UniProt ID:	P12849
Synonyms:	AI385716; Ribeta
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)