

Product datasheet for **MP211971**

Pofut2 Mouse qPCR Primer Pair (NM_030262)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	80294
Forward Sequence:	GAGACGAGTTTAGGAGCCAGCA
Reverse Sequence:	CCTAGTGCTGAGCCAGTTTGA
ACCN:	<u>NM_030262</u> , <u>NM_030262.1</u> , <u>NM_030262.3</u> , <u>BC147078</u> , <u>BC003494</u> , <u>BC018194</u> , <u>BC147077</u>
UniProt ID:	<u>Q8VHI3</u>
Synonyms:	231001IG23Rik; AI256847; BC003494; C2Iorf80; FUT13
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

