

## Product datasheet for **MP211999**

### PPP1R8 Mouse qPCR Primer Pair (NM\_146154)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	100336
Forward Sequence:	GAAGCCTCAGACATTGCCATCG
Reverse Sequence:	GGCAGTGTGAACTCTGTCAGG
ACCN:	<u><a href="#">BC025479</a></u> , <u><a href="#">NM_146154</a></u> , <u><a href="#">NM_146154.1</a></u> , <u><a href="#">NM_146154.2</a></u> , <u><a href="#">NM_146154.3</a></u> , <u><a href="#">BQ749177</a></u> , <u><a href="#">BY257668</a></u> , <u><a href="#">BY293277</a></u>
UniProt ID:	<u><a href="#">Q8R3G1</a></u>
Synonyms:	6330548N22Rik; AU044684; NIPP1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)