

Product datasheet for **MP212003**

Ppm1a Mouse qPCR Primer Pair (NM_008910)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	19042
Forward Sequence:	GGACAAGTACCTGGAGAGCAGA
Reverse Sequence:	GGGATGTTCTCACTGGCTAACG
ACCN:	BC008595 , NM_008910 , NM_008910.1 , NM_008910.2 , NM_008910.3 , BE861326 , BM932496 , BY739448
UniProt ID:	P49443
Synonyms:	2310003C21Rik; 2900017D14Rik; AI427932; AU017636; MMPa-2; MPPa-1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)