

Product datasheet for **MP212024**

PPP2R5C Mouse qPCR Primer Pair (NM_012023)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	26931
Forward Sequence:	GCCTCTTTCCGACAGTTAGCC
Reverse Sequence:	CAGAATCTCGCTGCGTTGTCAC
ACCN:	<u>BC003979</u> , <u>NM_012023</u> , <u>NM_012023.1</u> , <u>NM_012023.2</u> , <u>NM_012023.3</u> , <u>BC147179</u> , <u>BC147180</u> , <u>BQ175856</u> , <u>NM_012023.4</u>
UniProt ID:	<u>Q60996</u>
Synonyms:	2610043M05Rik; 2700063L20Rik; AI060890; AW545884; C85228; D12Bwg0916e; mKIAA0044
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)